MORCHER® Implants

Ophthalmic innovations since 1955

ANIRIDIA IOL

ANIRIDIA IOL

IRISMATCH® IOL	03
ANIRIDIA RINGS	04
SULCUS RING	05
PARTIAL ANIRIDIA RINGS	05
ANIRIDIA IMPLANTS	06
PARTIAL ANIRIDIA IMPLANTS	08
GENERAL INFORMATION	11

IRISMATCH® IOL

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The MORCHER® Irismatch® IOL is an Aniridia Implant made of PMMA, composed of a central 3.0 mm optic and a colored diaphragm with an iris pattern.

Selecting the color

The iris color is selected from the color samples, which are presented in our Irismatch® Color Case. The correct Iris color is chosen by holding the sample close to the patient's eye.

*The Irismatch® Color Case contains a set of a **45 artificial color match irises**. A case may be obtained from your distributor.*

Please note

*It is important to note in advance that the implanted device is never 100% identical to the natural color of the patient's iris. It is important to do the match-up in ambient light, not artificial light. The Irismatch® tends to appear lighter in the eye once implanted. This is especially true for lighter colors such as blue, for example. **Therefore, we recommend selecting one to two nuances darker!***


How to order

Each Irismatch® Implant is custom-made. The lead time for production is approximately eight weeks.

The following data must accompany your order:

- Reference of the Model (Type 30B)
- Selected color code (e.g. E010)
- Diopter (e.g. +22.0 D.)
- Patient (e.g. Mr. A. Miller)

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*In each case two of these aniridia rings are implanted in the capsular bag and dialled against one another to form a closed diaphragm. **After the implantation of the rings, the IOL is implanted.** The capsular bag has to be intact! The use of sufficient visco elastic makes the implantation easier.*

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ADVANTAGES


- *small incision*
- *additional stabilization of the capsular bag*
- *well centered*

TYPE 50C

Design	Rasch/Rosenthal
Total Diameter	10.0 mm
Diaphragm	Covers up to 12 cl. hours / 360°
Inner Diameter	6.0 mm
Incision	> 2.0 mm

TYPE 50E

Design	Rasch/Rosenthal/Zirm
Total Diameter	10.0 mm
Diaphragm	Covers up to 12 cl. hours / 360°
Inner Diameter	3.5 mm
Incision	> 3.25 mm

TYPE 50F

Design	Miller
Total Diameter	10.0 mm
Diaphragm	Covers up to 12 cl. hours / 360°
Inner Diameter	4.0 mm
Incision	> 3.0 mm

TYPE 96C NO


Sulcus Ring	
Design	Miller
Total Diameter	13.0 mm
Diaphragm (Aniridia Ring)	Covers up to 12 cl. hours / 360°
Diaphragm (Partial Aniridia Ring)	Covers up to 6 cl. hours / 180°
Inner Diameter	4.0 mm
Incision	> 4.5 mm

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Capsular tension ring with a diaphragm to dim a coloboma in the eye not larger than three hours. Is the coloboma three to six hours two rings can be implanted.

TYPE 96F

Design	Annen
Total Diameter	11.0 mm
Diaphragm	Covers up to 3 cl. hours / 90°
Inner Diameter	4.0 mm
Incision	> 3.5 mm

TYPE 96G

Design	Rasch / Knittel / Sold
Total Diameter	11.0 mm
Diaphragm	Covers up to 3 cl. hours / 90°
Inner Diameter	4.0 mm
Incision	> 3.5 mm

TYPE 96E


Design	Annen / Climenhaga
Total Diameter	11.0 mm
Diaphragm	Covers up to 3 cl. hours / 90°
Inner Diameter	5.5 mm
Incision	> 3.0 mm

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
Posterior Chamber IOL with diaphragm for idiopathic and traumatic aniridia.

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



TYPE 94A

Design	Sundmacher / Öhrström
Total Diameter	13.0 mm
Diaphragm	6.5 – 10.0 mm
Haptic Angulation	0°
Incision	> 10.0 mm
Diopter range	NONE (only available without optic)

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The ophthalmologist must choose each implant according to the state an size of the partial aniridia.

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In 1991, in cooperation with Prof. Sundmacher of the University of Düsseldorf, MORCHER® developed a special IOL with an integrated iris diaphragm to correct congenital and traumatic forms of aniridia.

Aniridia is a condition where the iris or part of the iris is missing and thus light falls on the retina unimpeded, blinding the individual. Aniridia can be congenital or result from injury. The implantation of an aniridia implant is to eliminate the glare caused by aniridia.

Our aniridia implants have been clinically tested since 1991 and were CE certified in the course of a product certification in 1994.

ADVANTAGES

- *By the implantation of aniridia implants a controlled amount of light enters the eye and, at the same time, the cataract may be treated.*
- *Product diversity*
- *Individual dioptre range*
- *Approved for more than 10 years*
- ***Custom-made designs***

PLEASE NOTE

When handling aniridia implants please note that black PMMA is more brittle than the transparent, pre-treated PMMA material. The crossing of the haptic to the diaphragm, which is especially brittle, is to be treated with the utmost care.

A-CONSTANTS

All numbers listed within the catalog are presented as guidelines only. We recommend you develop your own A-constants based on your experience.

Morcher GmbH

Kapuzinerweg 12
D-70374 Stuttgart

Telefon +49 (0) 711 / 95 320 - 0
Telefax +49 (0) 711 / 95 320 - 80

E-Mail info@morcher.de
Internet www.morcher.com