

**unique**  **edge**<sup>®</sup>  
MICROSURGICAL KNIVES



**Precision you can see. Quality you can feel.**



Clinicians know that minimizing risks and providing positive patient outcomes are the most important issues to consider when choosing the right surgical instruments. Unique Technologies offers a line of UniqueEdge<sup>®</sup> knives that consistently provide both the quality and precision necessary to achieve these goals.

### PRECISION YOU CAN SEE

- Our vacuum-purified double melt, hardened and tempered stainless steel technology produces a uniform strength and ductility that is critical to creating our precise cutting edge for smooth, uniform incisions
- Consistent blade dimensions
- Flawless cutting edge
- Superior wound architecture
- Minimal wound leakage

### QUALITY YOU CAN FEEL

- Effortless tissue penetration
- Excellent tactile response with minimal drag
- Consistent reproducible incisions

### PERFORMANCE YOU CAN TRUST


- Accurate, well controlled incisions
- Tightly sealing incisions
- Sutureless closure
- Rapid wound healing
- Reduced potential for infection

#### TEXTURED HANDLE

- Market proven design
- Enhanced instrument control
- Excellent tactile response

#### MATTE FINISHED BLADE

- Improves view of the operating field
- Eliminates reflective glare from operating room light sources
- Maximal visibility of the knife within the eye chamber





### PRECISION YOU CAN SEE

The microscope reveals the flawless, finished cutting edge that makes our microsurgical blades instruments of precision.

### BROAD CUSTOMIZATION CAPABILITIES

With the wide variety of blade dimensions and geometries available, the UniqueEdge® line is ready to meet your most demanding surgical requirements.

You can operate with complete confidence with a knife that flawlessly compliments your procedure, equipment, instruments, and surgical technique.



## Stab Knives



ITEM NO.	WIDTH	DESCRIPTION
52-1501	15°	Full Handle
52-2201	22.5°	Full Handle
52-3001	30°	Full Handle
52-4501	45°	Full Handle
52-1531	15° 3mm	Full Handle, restricted depth blade
52-1551	15° 5mm	Full Handle, restricted depth blade
52-1532	15° 3mm	Blade and tip, restricted depth blade
52-1552	15° 5mm	Blade and tip, restricted depth blade
52-1502	15°	Blade and tip
52-2202	22.5°	Blade and tip
52-2241	22.5° 4mm	Full Handle, restricted depth blade
52-3036	30° 3.5mm	Full Handle, restricted depth blade

## Clear Corneal Knives



ITEM NO.	WIDTH	DESCRIPTION
52-2563	2.5mm	Angled, full handle
52-2663	2.65mm	Angled, full handle
52-2763	2.75mm	Angled, full handle
52-2863	2.85mm	Angled, full handle
52-3063	3.0mm	Angled, full handle
52-3263	3.2mm	Angled, full handle

## Clear Corneal Knives- Side Cutting/ Parallel Sides



ITEM NO.	WIDTH	DESCRIPTION
52-2635	2.65mm	Angled, full handle
52-2735	2.75mm	Angled, full handle
52-2835	2.85mm	Angled, full handle
52-3035	3.0mm	Angled, full handle
52-3235	3.2mm	Angled, full handle

## Phaco Slit Knives



ITEM NO.	WIDTH	DESCRIPTION
52-2561	2.5mm	Angled, full handle, bevel up
52-2661	2.65mm	Angled, full handle, bevel up
52-2761	2.75mm	Angled, full handle, bevel up
52-2861	2.85mm	Angled, full handle, bevel up
52-2961	3.0mm	Angled, full handle, bevel up
52-3261	3.2mm	Angled, full handle, bevel up
52-3561	3.5mm	Angled, full handle, bevel up
52-5261	5.2mm	Angled, full handle, bevel up
52-2901	3.0mm	Straight, full handle, single bevel
52-3201	3.2mm	Straight, full handle, single bevel
52-2762	2.75mm	Angled, blade and tip, bevel up
52-2862	2.85mm	Angled, blade and tip, bevel up
52-2962	3.0mm	Angled, blade and tip, bevel up
52-3262	3.2mm	Angled, blade and tip, bevel up
52-2831	2.85mm	Angled, full handle, double bevel
52-2931	3.0mm	Angled, full handle, double bevel

## MicroPhaco Slit Knives



ITEM NO.	WIDTH	DESCRIPTION
52-1463	1.4mm	Angled, full handle, bevel up
52-1563	1.5mm	Angled, full handle, bevel up
52-1863	1.8mm	Angled, full handle, bevel up
52-2063	2.0mm	Angled, full handle, bevel up
52-2261	2.2mm	Angled, full handle, bevel up
52-2263	2.2mm	Angled, full handle, double bevel

Please call your UniqueEdge® distributor to inquire about knife designs or sizes not listed.

*Precision you can see. Quality you can feel.*

## Crescent Knives



ITEM NO.	WIDTH	DESCRIPTION
54-1060	1.25mm	Angled, full handle, bevel up
54-1010	2.0mm	Angled, full handle, bevel up
54-1000	2.0mm	Straight, full handle
54-5062	2.0mm	Angled, blade and tip, bevel up
54-5000	2.0mm	Straight, blade and tip
54-1013	2.3mm	Angled, full handle, bevel up

## Vitrectomy Knives



ITEM NO.	WIDTH	DESCRIPTION
51-1901	19 gauge	Straight, full handle
51-2001	20 gauge	Straight, full handle
51-2061	20 gauge	Angled, full handle, bevel up
51-2301	23 gauge	Straight, full handle
51-2361	23 gauge	Angled, full handle, bevel up
51-2501	25 gauge	Straight, full handle

## Trapezoid Knives



ITEM NO.	WIDTH	DESCRIPTION
55-1214	1.2mm–1.4mm	Angled, full handle
55-1416	1.4mm–1.6mm	Angled, full handle
55-1520	1.5mm–2.0mm	Angled, full handle
55-2724	2.75mm–3.2mm	Angled, full handle
55-2824	2.85mm–3.2mm	Angled, full handle

## Sideport Knives-MVR Design



ITEM NO.	WIDTH	DESCRIPTION
58-2008	0.8 mm	Angled, full handle
58-2010	1.0mm	Angled, full handle
58-2012	1.2mm	Angled, full handle
58-2014	1.4mm	Angled, full handle

## Implant Knives



ITEM NO.	WIDTH	DESCRIPTION
56-5361	3.2mm	Angled, full handle, bevel up
56-3561	3.5mm	Angled, full handle, bevel up
56-4061	4.0mm	Angled, full handle, bevel up
56-5261	5.2mm	Angled, full handle, bevel up
56-5561	5.5mm	Angled, full handle, bevel up
56-6261	6.2mm	Angled, full handle, bevel up

## Specialty Knives



ITEM NO.	DESCRIPTION
58-5700	Scleral, straight, blade only
58-6400	Round Tip, straight, blade only
58-6600	Lamellar, angled, blade only
58-6700	Mini, straight, blade only
58-6810	Needle blade, straight, full handle
58-6900	All Around, straight, blade only

## Disc/Spoon Knives



ITEM NO.	WIDTH	DESCRIPTION
54-2261	2.2mm	Angled, full handle, bevel up
54-3061	3.0mm	Angled, full handle, bevel up

## Myringotomy Blades



ITEM NO.	DESCRIPTION
58-7100	45°, downward cut, offset
58-7110	45°, upward cut, offset
58-7120	Spear Tip, narrow shaft, offset
58-7121	Lance Tip, narrow shaft, offset

To locate a UniqueEdge® Distributor or for ordering information, contact us at [www.unique-tech-inc.com](http://www.unique-tech-inc.com) or by Phone at 610-775-9191

### **Made in the USA**

- 100%
- Satisfaction guaranteed
- ISO and CE certification confirms consistency and reliability of manufacturing processes and product quality

### **Private label opportunities**

- We manufacture for some of the largest and most trusted names in the ophthalmic industry
- Contact us with your specific product needs

ISO 9001-2000 certified & ISO 13485-2003 certified

**CE** Mark Authorized under MDD 93/42/EEC

**unique**  **technologies inc.**

Mohnton PA USA

Office: 610-775-9191 Fax: 610-775-9192

 [www.unique-tech-inc.com](http://www.unique-tech-inc.com)