



**unique**  **edge**<sup>®</sup>  
MICROSURGICAL KNIVES

## Performance, Precision *and* Protection


**UniqueEdge offers a new line of guarded microsurgical knives that are simple to use, fully retractable and have a locking design that makes safety easy.**

- Uses the same handle that you've always valued
- Requires no change in your surgical technique
- Retains maximal feel and absolute control
- Protects your staff from accidental injuries
- Protects the blade from accidental damage
- Part of your OSHA Exposure Control Plan (ECP)



**Unique Edge® knives that provide quality and precision are now available with a simple, easy to use guard.**

To Use the Knife, Simply:



To Dispose of the Knife, Reverse Steps and Place in a Sharps Container


**Guarded Stab Knives**

ITEM NO.	WIDTH	DESCRIPTION
52-1501-G	15° Full Handle	
52-3001-G	30° Full Handle	



**Guarded Crescent Knives**

ITEM NO.	WIDTH	DESCRIPTION
54-1010-G	2.0mm	Angled, full handle, bevel up



**Guarded Phaco Slit Knives**

ITEM NO.	WIDTH	DESCRIPTION
52-2761-G	2.75mm	Angled, full handle, bevel up
52-2861-G	2.85mm	Angled, full handle, bevel up
52-2961-G	3.0mm	Angled, full handle, bevel up
52-3261-G	3.2mm	Angled, full handle, bevel up



**Guarded Sideport Knives - MVR Design**

ITEM NO.	WIDTH	DESCRIPTION
58-2010-G	1.0mm	Angled, full handle



**Guarded MicroPhaco Slit Knives**

ITEM NO.	WIDTH	DESCRIPTION
52-2263-G	2.2mm	Angled, full handle, double bevel



**TEXTURED HANDLE** .....

- Market proven design
- Enhanced instrument control
- Excellent tactile response

**MATTE FINISHED BLADE** .....

- Improves view of the operating field
- Eliminates reflective glare

**PROTECTIVE GUARD\*** .....

- Protects from sharps injuries
- No change in surgical technique
- Simple, easy to use



To locate a UniqueEdge® Distributor or for ordering information, contact us at [www.unique-tech-inc.com](http://www.unique-tech-inc.com) or by Phone at 610-775-9191

ISO 9001-2000 certified & ISO 13485-2003 certified  
**CE** Mark Authorized under MDD 93/42/EEC

\* patent pending

**unique**  **technologies inc.**

Mohnton PA USA Office: 610-775-9191 Fax: 610-775-9192